

SeaHARRE-5: lab E

«The Bodø method»

Einar Skarstad Egeland
Faculty of Biosciences and Aquaculture
Bodø University College
P.O.Box 1490
NO-8049 Bodø, Norway
einar.skarstad.egeland@hibo.no

Storage & extraction

- ❑ Samples stored for 10 months at -30°C in original packing
- ❑ Each freeze-dried filter torn to pieces and soaked in 1 ml cold (-30°C) 30 % methanol in acetone, back in freezer
- ❑ Extraction over night, filtered, analysed by HPLC at the same day (or weekend)

HPLC instrument

- Agilent 1100 HPLC with quaternary pump, cooled autosampler with enlarged injection loop, thermostatted column, diode array detector

HPLC method

- ❑ Injection: 50 μ l of 1 ml extract, without mixing
- ❑ Two columns: 2 x ACE C18 5 μ m, 250 x 4.6 mm with separate guard column
- ❑ Detection at 390, 420, 450 and 480 nm
- ❑ Calibration by injecting various volumes of the conc. pure standards

HPLC method – cont.

- ❑ Solvents: 1 M ammonium acetate, methanol, acetone, hexane
- ❑ Gradient over 130 minutes from 20 % ammonium acetate in methanol to 60 % hexane in acetone (to rinse the column after each sample)

Alternative method


- ❑ **Single column: Brownlee (Perkin Elmer) C18
5 μm , 220 x 4,6 mm with guard column**
- ❑ **Separation not as good as double ACE
column**
- ❑ **Analysed, but not summarised due to lack of
technician & lack of time**

Pigment separation


- ❑ Methods not checked for DV/MV-chl cs
- ❑ Not separating chl c2/MgDVP; DV/MV-chl bs; lutein/zeaxanthin
- ❑ Only alternative method separate chl c2/peridinin satisfactory
- ❑ Calibration at the best wavelength for each pigment (390/420/450/480 nm)

DHI-mix # 105

DAD1 A, Sig=390,30 Ref=900,100 (NEONV077-B501.D)
 DAD1 B, Sig=420,30 Ref=900,100 (NEONV077-B501.D)
 DAD1 C, Sig=450,30 Ref=900,100 (NEONV077-B501.D)
 DAD1 D, Sig=480,30 Ref=900,100 (NEONV077-B501.D)




Natural sample: AG09



□ (Sorry for wrong identities here)


Natural sample: K34



Chlorophylls a / b

- ❑ No DV-chl a observed in natural samples
- ❑ Both chl a and chl b separated into several peaks (epimer, allomer(s), ...)

Example: K34



Chlorophyll problems

- ❑ **Not observed in DHI standards after evaporation and re-dissolving**
- ❑ **Not observed in DHI mixes 104, 105, 106**
- ❑ **Observed in natural samples & Sigma standards**

Chl. problems – cont.

- ❑ Sigma standard found 98 % pure at HPL, 44 % in Bodø (if main peak chl a)
- ❑ Natural samples: Probably artefacts during storage/freeze-drying/extraction
- ❑ Not present in extracts at other labs?
- ❑ Or: Not shown in their HPLC methods?

